# Instructions

## for use for Electronic Handheld Magnifier *eMag 43*

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
Congratulations on your purchase of the eMag 43. This electronic handheld magnifier allows you to magnify printed text and objects from 2 to 14 times the object's original size. The eMag 43 features a 4.3-inch, wide-screen (16:9 aspect ratio) LCD display.

Information in this document is subject to change without notice. No part of this publication may be reproduced or transmitted in any form or any means electronic or mechanical, for any purpose, without the express written permission of A. SCHWEIZER GmbH.

Copyright © 2014 A. SCHWEIZER GmbH. All Rights Reserved.

LED RADIATION  
DO NOT STARE  
INTO BEAM  
CLASS 2 LED

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CH-8274 Tägerwilten

# Table of contents

<b>Safety and maintenance</b>	<b>6</b>		
Cleaning instructions	6		
<b>Quick start</b>	<b>7</b>		
<b>Overview</b>	<b>8</b>		
What's in the box?	10		
Batterie safety notes	10		
Inserting batteries	12		
Charging batteries	14		
LED descriptions	15		
<b>Using the eMag 43</b>	<b>16</b>		
Reading with the eMag 43	16		
		Magnifying an image	17
		Magnification level with handle collapsed	18
		Magnification level with handle extended	18
		Freezing and magnifying an image	19
		Saving and deleting an image	20
		Saving an image when the memory is full	21
		Turn LED Lights Off	21
		<b>Hardware description</b>	<b>22</b>
		<b>Carrying the eMag 43</b>	<b>23</b>
		<b>Service and support</b>	<b>24</b>
		Troubleshooting	24

# Safety and maintenance

To keep your eMag 43 magnifier in good condition, please observe the following guidelines:

Read the instructions in this user's guide before attempting to use the eMag 43.

Please do not attempt to repair the eMag 43 yourself or open its case. This product contains no user serviceable parts.

Keep the eMag 43 away from wet areas, rain, and other sources of liquids or environments with high humidity. Do not submerge the eMag 43 in water or other liquids.

Always operate your eMag 43 in a place where the temperature is between 10° and 35° Celsius.

Do not store the eMag 43 in places subject to temperature extremes such as in a parked car.

Do not attempt to use force when connecting the cables to the eMag 43. All appropriate connectors fit easily into the ports. Use only the power adapter that came with the eMag 43 to charge the NiMH batteries (nickel-metal hydride). Other adapters may damage the unit.

## Cleaning instructions

**CAUTION:** When cleaning the screen or camera window, turn off the unit, unplug the power cord, and apply a small amount of non-abrasive cleanser, such as the following, to a soft cloth (not a paper towel) and gently wipe: water, isopropyl alcohol or vinegar mixed with water at a concentration of no more than 10 % vinegar. Do not use any cleaner that contains the following agents: acetone, ethyl alcohol, ethyl acid, ammonia, or methyl chloride. Do not spray cleanser directly on the screen. Do not apply excessive pressure to the screen.


# Quick start

Note: Before using the eMag 43 electronic handheld magnifier, you must charge the unit. The batteries charge faster when the unit is off. Allow at least two hours for a complete charge before using the unit for the first time. Please observe the "Inserting and Charging the Batteries" chapter.


The eMag 43 is easy to use, even for those unfamiliar with assistive technology. Press the green/yellow **power**  button to turn the unit on. Hold or place the eMag 43 over an object and adjust the magnification to your preference using the yellow **magnification**  button. The magnified image appears on the unit's display.

For safety reasons, eMag 43 turns off automatically after 30 minutes if no button is pressed during that time. The following describes the location of the eMag 43 components and the LED operation.

# Overview

## What's in the box?

Before using the eMag 43, please take out and verify all the items that are in the box:

- eMag 43 electronic handheld magnifier
- AC power adapter with four international adapters
- Four AAA nickel-metal hydride (NiMH) rechargeable batteries
- Carrying case with shoulder strap
- Wrist strap
- eMag 43 User's Guide

## Charging batteries

Note: Before using the eMag 43 for the first time, you must charge the four AAA NiMH batteries.


Four AAA NiMH rechargeable batteries are included with the eMag 43. You can also use four micro type batteries (LR3, AAA, AM-4) to operate the unit. In this case, do not connect the power adapter.

Note: Do not mix NiMH batteries with alkaline batteries. Replace all batteries with four new ones of the same chemistry at the same time. Do not mix new batteries with old batteries. Failure to follow this precaution can lead to battery failure or leakage.  
Dispose of NiMH and alkaline batteries properly. They should not be disposed of with normal household waste. Check with your local government for disposal and recycling procedures.

CAUTION: Risk of malfunction if a NiMH or alkaline battery is replaced with an incorrect type.

We would like to inform you that our batteries/rechargeable batteries are marketed in accordance with the current law and refer you to the following points:

1. eMag 43 is operated by rechargeable batteries type AAA HR03. You can also use batteries type AAA LR03.
2. To prevent harmful effects to the environment and health, used batteries/rechargeable batteries must be disposed of properly at community recycling centres or returned to the merchant. Returning batteries is required by law and is completely free-of-charge for you.
3. Batteries/rechargeable batteries should only be returned in a discharged state.
4. Valuable raw materials such as zinc, iron and nickel can be reclaimed by recycling used batteries/rechargeable batteries.
5. The batteries/rechargeable batteries can be easily removed from our devices for disposal (see User's Guide).
6. The symbol of a "crossed-out trash container" indicates that it is illegal to dispose of used batteries with regular domestic waste.


The additional specification of the chemical sign for lead (Pb), cadmium (Cd) and mercury (Hg) indicates that these materials are contained in the batteries/rechargeable batteries in accordance with the legally-specified limits.


## Inserting batteries

To insert the batteries, do the following:


1. Hold the eMag 43 in one hand. With the other hand, rotate and swing the handle until it is fully extended as shown in the following illustration.

2. Turn the eMag 43 so that the LCD display is facing away from you, then slide the battery cover in the direction of the arrow shown in the following illustration.

Note: Observe the correct polarity of the batteries when inserting. Insert the batteries as shown in the following illustration.


3. Close the battery cover again.

## Charging batteries

To charge NiMH batteries, do the following:

1. Plug the small end of the power adapter into the power port. The port is located on the left top-edge of the unit next to the LED as shown in the following illustration.

2. Plug the other end of the power adapter into a standard AC outlet. The battery indicator LED blinks green to show that the batteries are charging.
3. Once the batteries are charged, the LED stops blinking and glows a steady green. Unplug the power adapter from the eMag 43. The battery indicator LED turns off.

## LED Descriptions

Red – Flashing	Battery low
Green – Flashing	Battery charging - power adapter connected
Green – Solid	Battery charged - power adapter connected
Red – Solid	Error – unplug power adapter, wait a few seconds and plug in again
No LED	Battery charge ok, power adapter not connected or the unit is off

Your eMag 43 is now ready for up to two hours of continuous use. When the batteries are low, the battery indicator LED flashes red. Connect the **power adapter** and plug it in to recharge the batteries.

Note: Please note that after 30 minutes the eMag 43 switches to energy-saving mode (screen goes blank) if no key is pressed during this time.  
Press the On/Off Button to reactivate the eMag 43.

# Using the eMag 43

After charging the batteries, you are ready to begin reading with the eMag 43. With the unit's powerful magnification features, you can read or view letters, menus, maps, photos, labels, bus and train schedules, television listings, receipts, ingredients, and much more.

GB

## Reading with the eMag 43

To read by using the eMag 43, do the following:

1. Make sure that the unit is charged.
2. Press the green/yellow **power**  button to turn on the eMag 43. After a few seconds, the LCD displays the current view as seen through the camera window.
3. If you are reading a book, letter, or other printed item, place it on a flat surface, such as a table or desk.
4. You can then hold the eMag 43 by its handle the same way as you would hold a normal magnifying glass, or you can fold the handle under the unit and place the eMag 43 directly over the text or object that you want to read or view. Anything under the unit is magnified and shown on the screen.
5. Adjust the zoom using the yellow **magnification**  button located on the left side next to the LCD.
6. Select a colour mode by pressing the blue **colour select**  button located on the left side of the screen. Five different video modes are available: Normal (full colour), White on Black (enhanced positive), Black on White (enhanced negative), Yellow on Blue and Yellow on Black.
7. Move or slide the unit to the right to read to the end of the line. Slide the unit first to the left and then down to read the next line of text.

For safety reasons, eMag 43 turns off automatically after 30 minutes if no button is pressed during that time.

## Magnifying an Image

The eMag 43 has three digital zoom settings. Depending on the position and placement of the unit, it supports screen magnification levels from **2x** to **14x**.

1. To increase screen magnification, press the yellow **magnification +** button.
2. To step up to the next level of magnification, press the yellow **magnification +** button again.
3. Pressing the yellow **magnification +** button a third time returns the eMag 43 to its initial magnification level.

For more information about the different magnification levels, please refer to the following section.

## Magnification level with handle collapsed


When the eMag 43 is placed on a document, the starting magnification is **5x**.

By pressing the yellow **magnification +** button repeatedly, the magnification increases to **7.5x**, to **10x**, and then back to **5X**.

## Magnification level with handle extended


When the eMag 43 is held over a document (as shown in the illustration), the starting magnification is **2x**.

Magnification level **2x** while holding the eMag 43 approximately 7.6 cm above item.

By pressing the yellow **magnification +** button repeatedly, the magnification increases to **3x**, to **4x**, and then back to **2x**.

When the eMag 43 is in this low position (approx. 1.3 cm) above a document, the starting magnification is **7x**.

By pressing the yellow **magnification +** button repeatedly, the magnification increases to **10.5x**, to **14x**, and then back to **7x**.

### Freezing and magnifying an Image

Freezing an image can be useful for reading labels on items on a store shelf, for example. By holding the eMag 43 close to the item and capturing the image, you can then hold the unit close enough for you to read.

To freeze the current image on the LCD screen, do the following:

1. Press the red **freeze frame**  button.
2. If you wish, you can change the magnification of the picture taken by pressing the yellow **magnification +** button.
3. Press the red **freeze frame**  button again to turn off the freeze frame mode.

## Saving and deleting an image

Switch on the eMag 43 (press the green/yellow button ). Then press the red button  to take the desired image. The image on the display can be saved as follows

- Press  and **HOLD this red button for 6 seconds** to enter Storage mode.
- The check box  in the lower left corner will be automatically checked. This indicates that the image will be stored when you exit Storage mode.
- The Forwards symbol is now displayed in the bottom right corner  Press the red button  to advance to the next image. Press the button  repeatedly to advance through all the images saved in the **eMag 43**
- To delete an image, press the blue button  to clear the check box The image is deleted when you exit Storage mode.
- Press and HOLD the red button  for **3 seconds** to leave Storage mode.
- To view stored images press and HOLD the red  button for **6 seconds**. You will then enter Storage mode You can then advance through the images as previously described in this section.

## Saving an image when the memory is full

The eMag 43 can store a maximum of 15 images. If you try to store another image, the "Memory full"  icon is displayed in the bottom left corner. Before you can save another image you must delete one or more images from the memory:


- Press the red button  to advance to the desired image and then press the blue button  to clear the check box.
- Press and HOLD the red button  for **3 seconds** to leave Storage mode. The new image is saved, and all unchecked images are deleted.

## Turn LED Lights Off

To avoid glare when viewing reflective surfaces you can turn off the LED lights.


- Press and HOLD  the green/yellow button and then press the red button . The LED lights turn off.
- To turn on the lights again, press and HOLD the green/yellow button  and then press the red button .

# Hardware Description

eMag 43 front and top view

1. Power button (green)
2. Colour select button (blue)
3. Magnification button (yellow)
4. Freeze frame button (red)
5. LCD screen
6. Power port
7. Battery indicator LED
8. Handle and battery compartment
9. Lights
10. Camera window
11. Wrist strap post


eMag 43 rear view

# Carrying the eMag 43

The eMag 43 has a wrist strap that you can use to carry it comfortably and safely. To attach the strap, do the following:

1. Push one end of the strap through the hole in the wrist strap post that is located on the rear of the eMag 43.
2. Loop one end of the wrist strap through the other end of the wrist strap and pull the strap tight.

When using the eMag 43, you can now place your hand in the wrist strap. This helps protect the eMag 43 from falling if you let it go by accident.

# Service and Support

Note: This device has no user-serviceable components. Any unauthorised attempt to service or replace internal components will void the product warranty.

GB

Please contact your merchant if you need technical assistance for your eMag 43.

## Troubleshooting

Problem	Solution
The screen is blank.	Ensure that the power is on. Verify that the item you want to magnify is positioned correctly under the camera window. Reduce the magnification level. Charge the NiMH batteries. If using AAA alkaline batteries, replace the alkaline batteries. (Do not mix NiMH batteries and alkaline batteries.)
The unit does not turn on.	Charge the NiMH batteries. If using AAA alkaline batteries, replace the alkaline batteries. (Do not mix NiMH batteries and alkaline batteries.)
The image on the screen is dark or distorted.	Make sure that the eMag 43 is positioned flat on the page you are reading.

The battery indicator LED is flashing red.	Charge the NiMH batteries, or, if using AAA alkaline batteries, replace the alkaline batteries. (Do not mix NiMH and alkaline batteries.)
The battery LED is solid red while the power adapter is connected to the eMag 43.	A battery error has occurred. Unplug the power adapter from the eMag 43 and plug it in again. Charge the NiMH batteries for at least 2 hours. If using AAA alkaline batteries, replace the alkaline batteries. If the problem continues, disconnect the power adapter and call your eMag 43 merchant.
The battery indicator LED does not flash green, even though the power adapter is connected and plugged into a wall socket.	Ensure that the small end of the power adapter is properly plugged into the unit's power port. Ensure that the outlet has power and is not malfunctioning. If the wall outlet is connected to a wall switch, make sure that the switch is on.
The NiMH batteries are not charging.	Ensure that the small end of the power adapter is properly plugged into the unit's power port. Ensure that the batteries are inserted properly. Observe the correct polarity. Ensure that the outlet has power and is not malfunctioning. If the wall outlet is connected to a wall switch, make sure that the switch is on.
There are smudges or blemishes on the screen.	Clean the LCD display and the camera window.
The text is too small to read.	Press the yellow <b>Magnification button</b> to increase the magnification level.
The text is too large to read.	Press the yellow <b>Magnification button</b> to cycle through to a lower magnification level.

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Copyright © 2014 A. SCHWEIZER GmbH. Tous droits réservés.
Copyright © 2011 A. SCHWEIZER GmbH. Alle rechten voorbehouden.